


# C86

## Patient Monitor



• 15 inch LED touch screen.

C86 patient monitor is designed to meet your every second care of patients in clinical, configuring 12.1" LED touch screen, fixed handle, various mounting solutions as well as handwriting pen, it is therefore your optimal choice for acute care. In case of different clinical environment such as in ICU, C86 provides IPX1 waterproof protection to satisfying strict environment requirements.



Fixed handle, more compact with small weight, easy to carry




Aesthetically pleasing new interface design



USB, VGA, network and multifunctional interface



Wall mount, rolling stand



Large capacity of Lithium battery support long time working without power supply

## Configuration

### Standard configuration

- ECG
- PR
- SpO<sub>2</sub>
- NIBP
- RESP


### Option configuration

- Dual-IBP
- BIS
- AG
- ICG
- Thermal Printer
- Nellcor/ Masimo SpO<sub>2</sub>
- EtCO<sub>2</sub>
- C.O.
- TEMP

## Recommended configuration

	operating room	ICU	CCU	general ward
12-lead ECG	✓	✓	✓	×
MASIMO SpO <sub>2</sub>	✓	✓	✓	×
Comen SpO <sub>2</sub>	✓	✓	✓	✓
Side-stream EtCO <sub>2</sub>	✓	✓	✓	×
Mainstream EtCO <sub>2</sub>	✓	✓	×	×
BIS	✓	✓	✓	×
C.O.	✓	✓	✓	×
IBP	✓	✓	×	×
AG	✓	×	×	×
ICG	✓	✓	✓	×

## COMEN



# C86

## Patient Monitor




### Shenzhen Comen Medical Instruments Co., Ltd.

Add: No.2 of FIYTA Timepiece Building, Nanhuan Avenue, Gongming Sub-district, Guangming New District, Shenzhen, P.R.China  
 Tel:+86-755-2640 8879 2641 9446 Fax:+86-755-2643 1232 Website: en.comen.com E-mail: info@szcomen.com

With leading ECG technology, anti-motion & weak perfusion SPO<sub>2</sub> technology as well as accurate NIBP measurement technology and cooperation with world leading medical technique providers such as Masimo, Covidien, Respironics, Medis, C86 is designed to optimize performances by configuring EtCO<sub>2</sub>, AG, BIS and noninvasive hemodynamic monitoring into one, helping you care even the most critical patients with professional assistance.


## ECG

- 3/5/12-lead ECG measurement technology, leads automatic identification
- Intelligent leads off detection and automatically leads selection guarantee uninterrupted monitoring
- ECG ensures intensive monitoring for a particular waveform
- CMRR  $\geq$  105dB, outstanding ECG anti-interference capability
- Support arrhythmia analysis & ST segment analysis




## IBP

2-channel IBP with SIMILAC accessories optional, Support monitor ABP, PAP, CVP, LAP, RAP, ICP etc




## NIBP

AcuTec™ NIBP technology, high accuracy for hypertension monitoring. The initial inflatable pressure can be selected to improve the accuracy of measurement and the comfort of patients.



## Anesthetic Gas

Collaborates with MASIMO, adopts the advanced anesthetic gas module for monitoring 8 types of gas (O<sub>2</sub>, Co<sub>2</sub>, N<sub>2</sub>O, ENF, ISO, DES, SEV, HAL). Automatic identification of the anesthetic gas, short time for warm-up, long service life and supports the MAC value (minimum alveolar concentration).



## EtCO<sub>2</sub>

- Collaborates with US RESPIRONICS, MASIMO, Plug and Play EtCO<sub>2</sub> monitoring.
- Use CAPNOSTAT 5 / IRMA mainstream sensor for optimal performance in monitoring intubated patient.
- Small, durable and lightweight mainstream sensor provides accurate and reliable monitoring for all intubated patients from neonates to adults.
- No calibration required.
- Use LoFlo / ISA sidestream sensor for monitoring non-intubated patient.
- Flexible, compact CO<sub>2</sub> sensor provides consistent and reliable monitoring of adult, pediatric and neonatal patients.
- Sample rate  $\leq$  50ml/min(micro-stream).




## BIS (Depth of Anesthesia) module

The BIS module has been designed to be used in the monitoring of the level of consciousness of a person during the application of general anaesthesia or in intensive care. This is accomplished by registering the electroencephalographic signal (EEG) by means of surface electrodes which is then analyzed by a digital process. As a result of the applied calculation, an index "BIS" is obtained, which serves as guidance to the experts who use it to determine the level of consciousness of the patient during surgery.

# Critical time Seize every second to save life

## Non-invasive Hemodynamic

- Collaborates with MEDIS, impedance cardiography for non-invasive continuous hemodynamic monitoring.
- Micro-signal transmit through disposable electrode.
- Blood volume and Blood Flow Velocity varies with heartbeat, DISQ® technology processes impedance signal variation.
- Variation of impedance applies to non-invasive Z MARC™ algorithm for acquiring SV, CO, SVR, Contractility and TFC etc.



## Intelligent Alarm


I-KLOK® intelligent alarm management, auto-identification of alarm level. Self-adjust proper alarm time to reduce false alarms.

## C.O. (invasive cardiac output) module

C86 is involved itself in invasive cardiac output technique, but C.O. measurement is conducted with conventional thermo dilution invasive cardiac output and other hemodynamic parameters. The monitor can measure "blood temperature", "calculating cardiac output", "calculating hemodynamics". The cardiac output is measured with floating catheter led from vein to pulmonary artery followed by injecting a certain amount of ice water at 0°C (injecta) such that the blood temperature will be varied after the injecta and blood output from the heart are mixed together thereby achieving cardiac output by measuring blood temperature variation before and after injected in accordance with the principle of heat balance.

## Masimo SpO<sub>2</sub>

Performance Claim	MasimoSET Pulse Oximeter
<b>SpO<sub>2</sub> Accuracy (70-100%)</b>	
Adult/Pediatric (No Motion)	±2 digits
Adult/Pediatric (Motion)	±3 digits
<b>Perfusion Index Range</b>	0.02% - 20%
<b>Accuracy in Low Perfusion</b>	Adult ±2 Neo ±3 digits
<b>Forehead Sensor</b>	TF-1 ±2 digits
<b>Ear Sensor</b>	TC-1 ±3.5 digits
<b>Fragile Skin non-adhesive (No Motion)</b>	SoftTouch ±3 digits
<b>Fragile Skin non-adhesive (Motion)</b>	SoftTouch ±3 digits
<b>SpO<sub>2</sub> Accuracy (60-80%)</b>	
Adult/Pediatric (No Motion)	Not Currently Claimed <sup>2</sup>
<b>Forehead Sensor</b>	Not Available
<b>Pulse Rate Accuracy (70%-100%)</b>	
<b>Pulse Rate (No Motion)</b>	25 - 240 bpm ±3 digits
<b>Pulse Rate (Motion)</b>	25 - 240 bpm ±5 digits
<b>Pulse Rate - Low Perfusion</b>	25 - 240 bpm ±3 digits



- Support wire & wireless central monitoring system.