



Aeon7200A

Full Configuration

Anesthesia Machine

CE 0123

AEOMED
Reliable Quality Thoughtful Service

- With CE certificate, meets EU clinical safety requirement
- Large 8.4" LCD color touchscreen, monitoring of waveforms and parameters in different colors
- IPPV, Pressure, and SIMV modes can easily meet clinical requirements
- 6-tubes of flowmeters for 3-gas sources and suitable for various operation requirements
- Spirometry loops, for easier interpretation of ventilation effectiveness
- Optional multi-gas monitoring provides more reference for anesthesia
- Optional AGSS provides more protection

Technical Specifications

Physical & Environmental Specifications		Setting: Touchscreen and Navigator knob	
• Dimensions(H×W×D):	1356×845×620 mm	• Tidal volume:	20~1500 ml
• Weight:	110 kg (basic unit)	• Ventilation frequency:	2~100 bpm
• Storage temperature:	-20~55°C	• I:E ratio:	4:1~1:8
• Working temperature:	10~40°C	• PEEP:	0~30 cmH ₂ O (optional)
• Storage humidity:	≤93%	• Inspiratory plateau:	OFF, 5%~60%
Electrical Specifications		• Flow trigger:	1~15 L/min
• Mains:	AC 100~240 V, 50/60 Hz	• Pressure trigger:	-20~0 cmH ₂ O
• Battery:	DC 24 V, minimum 120min	• Pressure range:	5~70 cmH ₂ O
• Mains outlet:	3 (1.5A individual)	• Monitoring display:	Pressure: P _{peak} , P _{plat} , P _{mean} , PEEP
Pneumatic Specifications			Volumes: MV, Vt
• Gas supply:	0.28~0.6 MPa		Freq.
• Flowmeters:	O ₂ : 0~1 L/min; 1~10 L/min		Compliance
	N ₂ O: 0~1 L/min; 1~10 L/min		Gas: O ₂ , CO ₂ , Agent (all optional)
	Air: 0~1 L/min; 1~12 L/min	• Graph display:	P-t, F-t waveforms, P-V, F-V loops
• Gas system:	Oxygen supply failure protection	• Alarm: Audible and Visual alarms	
• O ₂ flush:	25~75 L/min		High/Low Airway Pressure
• Working mode:	Closed, Semi-closed, Semi-open		High/Low Minute Volume
• Driven mode:	Pneumatically driven and electronically controlled		High/Low FiO ₂ (when FiO ₂ function is enabled)
• Safety valve:	≤10 kPa		Power failure
• Operating mode:	One key Man. / Vent. switch		O ₂ Supply failure, Apnea,
• Cylinder yokes:	Optional	• Alarm Silence:	≤120 s
Ventilator Specifications		Vaporizer	
• Patient type:	Adult, Pediatric	• Supports 2 vaporizers (Selectatec® with interlock, optional standby vaporizer parking holder)	
• Ventilation modes: IPPV, PLV, SIMV, Spont (with apnea backup), Manual, Standby		• Agent type: Halothane, Enflurane, Isoflurane, Sevoflurane	
		• Filling type: Pour-fill, Key-fill, Quik-fil®	

Remark: Above configurations include standard and option. Please check price with your Aeonmed sales representative.

AEOMED
Reliable Quality Thoughtful Service

HQ: Building 9, No.26 Outer Ring West Road, Fengtai District, Beijing 100070, China
Science Park: No. 10, Chaobai Street, Yanjiao Development Zone, Sanhe City, Hebei Province
065201, China
TEL: +86-10-5841 1198 FAX: +86-10-6371 8989
Http: //www.aeonmed.com E-mail: marketing@aeonmed.com

The information contained in this leaflet is correct at the date of publication. The policy of Beijing Aeonmed Co., Ltd. is one of the continued improvement to its products. Because of this policy, Aeonmed reserves the right to make any change, which may affect the information in this leaflet without giving prior notices.

AAM7200A-1902