



Feature

- Colorful 15 inch TFT screen, multi-parameters, easy to observe.
- New ultra-thin wall-mounted style, space-saving ease of placement and operation.
- Touch screen operation, more efficient and convenient; and support the remote control operation.
- Interactive man-machine interface, multi-level menu, ample setup and displays the abundance intuitionistic contents.
- Have a variety of styles of multi-lead ECG, Big font, RESP breathing, oxygenation, short trend display interface.
- ECG、NIBP、SP02、RESP、TEMP、HR/PR
- EtCO2 (optional)

Functions:

- With user information management function, to input patient information and identity management.
- Comply with the YY709 alarm system, eye-catching sub-column independent warning lights can prompt to three sound and light alarm.
- With alarm, multi-level pulse volume adjustable function; the parameter display colors can be adjusted based on user choice.
- Freezing and releasing waveform function, analysis of ST-segment. With 192-hour trend graph, 1000 sets of NIBP, list, arrhythmia analysis and other functions.
- Dual-temperature measurement, can simultaneously measure the temperature and temperature difference between the display surface and the body cavity.
- Support the drug concentration in titration calculations and event functions customized and record.
- Have a unique infrared remote control function to facilitate long-distance rounds operation.
- With pacing analysis, can be used simultaneously with pacemakers, defibrillators.



ChenWei 晨伟

Configuration:

- Support DC power input, and adapt to a variety of environmental conditions.
- External replaceable backup lithium battery, the battery can be used more than two hours.
- Support external EtCO2 interface, can be connected to the mainstream and side stream CO2 module, to facilitate the clinical needs.
- 32 TCP/IP network connection, composed of central information system.
- External recorder can print waveform, data and other information in real time.
- External mobile transmitter module can be connected to central monitoring system.
- Desktop mounting bracket, two channel TEMP, tow channel IBP, wireless networking module, GPRS transmitter module.
- Can be fixed on wall, desktop, carts, ambulance, and others.

Nanjing Chenwei Medical Equipment Co., Ltd
Formerly Shenzhen Chenwei Electronic Co., Ltd
93 Shengtai West Road, Jiangning Development Zone, Nanjing China 211100
Tel: +86 25 6603 5699
Fax: +86 25 5791 8628
Email: sales@chenwei-med.com
Find us: www.chenwei-med.com

Authorized Representative

Note: specification subject to change without prior notice