90105201

GENERAL SPECIFICATIONS

- Laminate top is resistant to water and heat. Useful drawer.
- Long lived hinge and rail system. Strong drawer and cabined handle.
- Adjustable bed table according to height of bed and hidden beside of cabinet if necessary.
- Cavity designed for multi usage.
- Epoxy painted steel construction is resistant to corrosion.
- With four castors (two of them with brake) can move freely.

TECHNICAL SPECIFICATIONS

Width	:	50 cm
Depth	:	40 cm
 Height 	:	80,5 cm
Dining table height	:	85 - 110 cm
 Weight 	:	~ 20 kg
Volume	:	0,19 m ³

* RS : Melamine color options page 81. * RL : Laminate color options page 80.



BED SIDE CABINET

DESCRIPTION OF OPTION CODES KL1 LAMINATE COATED CHIPBOARD TOP MELAMINE COATED CHIPBOARD TOP KL4

Must select only one of them.



DOLSAN



BED SIDE CABINET

GENERAL SPECIFICATIONS

90105301

- Wood body is durable against to water and heat.
- Useful drawer.
- Long lived hinge and rail system.
- Strong drawer and cabined handle.
- Laminated drawer and door are durable against to water and heat.
- With four castors (two of them with brake) can move freely.

TECHNICAL SPECIFICATIONS

 Width 	: 40 cm
 Depth 	: 40 cm
 Height 	: 78,5 cm
 Weight 	: ~ 17 kg
 Volume 	: 0,13 m ³

* RS : Melamine color options page 81. * RL : Laminate color options page 80.



T8

KL1

KL4

DESCRIPTION OF OPTION CODES T8 PLASTIC CASTORS TWO WITH BRAKE T9 WITH ALUMINIUM FIXED LEG



