

90105101

GENERAL SPECIFICATIONS

- Laminate or ABS top is resistant to water and heat.
- Useful drawer.
- Long lived hinge and rail system.
- Strong drawer and cabined handle.
- A shelf designed for multi usage in the cabinet.
- Epoxy painted steel construction is resistant to corrosion.
- With four castors (two of them with brake) can move freely.

TECHNICAL SPECIFICATIONS

■ Width	: 40 cm
■ Depth	: 40 cm
■ Height	: 79,5 cm
■ Weight	: ~ 12 kg
■ Volume	: 0,14 m ³


* RS : Melamine color options page 81.

* RL : Laminate color options page 80.

BED SIDE CABINET

OPTIONS

GROUP 1 / CASTOR OPTIONS



GROUP 2 / TABLE TOP OPTIONS



DESCRIPTION OF OPTION CODES

- T8** PLASTIC CASTORS TWO WITH BRAKE (50 MM)
- T9** WITH ALUMINIUM FIXED LEG
- KL1** LAMINATE COATED CHIPBOARD TOP
- KL2** ABS TOP
- KL4** MELAMINE COATED CHIPBOARD TOP

Must select only one of them.



90105102

GENERAL SPECIFICATIONS

- Laminate or ABS top is resistant to water and heat.
- Useful drawer.
- Long lived hinge and rail system.
- Strong drawer and cabined handle.
- A shelf designed for multi usage.
- Epoxy painted steel construction is resistant to corrosion.
- With four castors (two of them with brake) can move freely.

TECHNICAL SPECIFICATIONS

■ Width	: 40 cm
■ Depth	: 40 cm
■ Height	: 78,5 cm
■ Weight	: ~ 12 kg
■ Volume	: 0,13 m ³

* RS : Melamine color options page 81.

* RL : Laminate color options page 80.

OPTIONS

GROUP 1 / CASTOR OPTIONS



GROUP 2 / TABLE TOP OPTIONS



DESCRIPTION OF OPTION CODES

- T8** PLASTIC CASTORS TWO WITH BRAKE (50 MM)
- T9** WITH ALUMINIUM FIXED LEG
- KL1** LAMINATE COATED CHIPBOARD TOP
- KL2** ABS TOP
- KL4** MELAMINE COATED CHIPBOARD TOP

Must select only one of them.

