

90102101 CONFIDENCE 1

HOSPITAL BED FOR CHILD

GENERAL SPECIFICATIONS

- Easily cleanable metal bed platform designed for the ventilation of the bed.
- Backrest adjusted by mechanical system.
- Ergonomic, durable and easily cleanable rails system that can be adjustable at both sides of the bed.
- Shock resistance plastic bumper wheels.
- Metal components dyed with anti-corrosion electrostatic powder paint after being treated with phosphate.
- Cabinet, which is metal bodied, is designed to multi usage on legrest of child bed. (Optional)
- Optional laminate color choice.^(RLX)

TECHNICAL SPECIFICATIONS


■ Width	:	77 cm
■ Length	:	156 cm
■ Height	:	50 cm
■ Weight (Without Accessory)	:	~ 53 kg
■ Accompaniment Weight	:	~ 13 kg
■ Weight (Cabinet)	:	~ 19,5 kg
■ Volume (Without Accessory)	:	0,17 m ³
■ Volume (Cabinet)	:	0,14 m ³

OPTIONS

GROUP 1 / CASTOR OPTIONS



GROUP 2 / CABINET OPTIONS



GROUP 3 / ACCOMPANIMENT OPTIONS



GROUP 4 / IV STAND OPTIONS



DESCRIPTION OF OPTION CODES

- T0** WITHOUT CASTORS
 - T4** METAL CASTOR TWO WITH BRAKE (125 MM)
 - D0** WITHOUT CABINET
 - D2** WITH CABINET
 - R0** WITHOUT ACCOMPANIMENT BED
 - R1** WITH ACCOMPANIMENT BED
 - W0** WITHOUT OTHER ACCESSORIES
 - W8** IV STAND : ONE HEIGHT ADJUSTMENT SYSTEM
- Must select only one of them.



* RL : Laminate color options page 80.