



CWM-301C Anesthesia system

Feature

- Integrated breathing circuit with APL valve, excellent tightness.
- Replacing drive bellow is NOT required when operating adult/pediatric anesthesia.
- Applied with 1 high-precision vaporizers. (Selectatec).
- 10.1 inch touch screen with no dead space.
- Real-time detection for loops, Paw, BPM, VT and compliance.
- Comprehensive monitoring capabilities include ventilation parameters monitoring, preset parameters, alarm indication and breathing waveform monitoring.
- Head light.
- Altitude compensation function.
- ACGO available to connect with T/Bain circuit for open operation(Optional).
- Yoke system with up to 3 back up cylinders (A type with 11L size) (Optional).




Technique Specifications	
Physical Specifications	
Dimensions:	70*85*140cm
Weight:	65kg
Casters:	Front wheels with lock
Screen:	10.1 inch LCD
Working Specifications	
Power	AC 100V-240v, 50Hz±2%
Pipeline Supply	
Gas configuration	O ₂ , N ₂ O, Air(Optional)
Pipeline input range	280-600kPa
Battery Power	
Run-time	≥180 mins
Pre-setting Functions	
Language	English
Calibration	Automatic
Applications	
Patient range	Audlt/ Pediatric
Ventilator Specifications	
Modes of Ventilation	
A/C	Assist / control Ventilation
VCV	Volume Controlled Ventilation
SIMV	Synchronized Intermittent Mandatory Ventilation
MANUAL	Manual
PCV	Pressure Control Ventilation
STAND-BY	STAND-BY
Ventilator Parameter	
Working type	Electronically controlled, pneumatic driven
Tidal Volume range	20-1500ml
Pressure range	5-60 cmH ₂ O
Rispiration rate range	1-100 bpm (1-40bpm under SIMV)
I:E range	4:1-1:8
Ptrigger range	-20-20 cmH ₂ O (based on PEEP)
Minute volume range	>18L/Min
SIGH range	0-5/100
System Standard	
Hypoxic guard system	N ₂ O cut-off valve, O ₂ concentration >25%
Safety valve	<12.5kPa
Fresh gas compensation	25-75L/min
Volume of CO ₂ absorber	2 L
Flow meter	4/5(Optional) tubes, Cascade 0.1-15L/min O ₂ 0.1-15L/min N ₂ O, 0.1-15L/min air
Vaporizer	1/2(Optional) Selectetac, maximum 2 mounts(Isoflurane/Enflurane/Sevoflurane/Halothane)
Monitoring	
Display graphics	Waveforms of P-T, F-T, V-T, loops of P-V,V-F,F-P
Types	VT, MV, BPM, Paw, Lung compliance, Inspiratory Platform,FiO ₂
Alarm Limits	
Audio/Visual	No tidal volume, MV,Paw, FiO ₂ limitation O ₂ , Air failure, AC power failure, Battery low
PEEP	
Type	Integrated, Electronic control
Range	OFF, 3-30cmH ₂ O
Optional	
1. AGSS	
2. Suction Unit	
3.Communication Interface	



Nanjing Chenwei Medical Equipment Co., Ltd
Formerly Shenzhen Chenwei Electronic Co., Ltd
93 Shengtai West Road, Jiangning Development Zone, Nanjing China 211100
Tel: +86 25 6603 5699
Fax: +86 25 5791 8628
Email: sales@chenwei-med.com
Find us: www.chenwei-med.com

Authorised Representative



Note: specification subject to change without prior notice