

Creation With Hearts

DF-200E

Defibrillator Monitor



Two Defibrillation Mode: AED/Manual

Advanced Biphasic Defibrillation Technology

promed

DF-200E



Display

Screen Type: High-resolution display
Screen Size: 7 inches(17.7 cm) diagonally
Sweep Speed: 25mm/sec
Information: Heart Rate, Lead/Pads, Alarm On/Off, SpO2, AED Functions and Prompts, Alarm Selection and Limits, Delivered Energy

Defibrillator

Waveform: Biphasic
Energy Display: Monitor display indicates both selected and delivered energy
Charge Prompt Type: Voice and visual prompts
Electrode Impedance Measurement Range: 0-250 ohms

ECG Monitoring

Patient Connection: 3-lead ECG cable, or 5-lead ECG cable, paddles
Lead Selection: Displayed on monitor, paddles, I, II, III, AVR, AVL, AVF, V
ECG Size: 0.25, 0.5, 1, 1.5, 2, 4 cm/mV display on monitor.
Heart Rate: 20-300BPM
Heart Rate Alarm: On/Off displayed on monitor, user-selectable
Smart Alarms: Beeper/voice prompts indicate shockable rhythm

Battery

Type: Rechargeable, Ni-MH battery, 12V
Operating Time: For a new, fully charged battery 60 defibrillator discharges at maximum energy, or 3 hours minimum of continuous ECG monitoring. Additional parameters will effect operating time with different functions.

NIBP Module (Optional)

Measurement Unit: mmHg/kPa
Measurement Range: Adult 10-270mmHg
Pediatric: 10-200mmHg; Neonate: 10-135mmHg

Sales Add: Room 705, Jinfulai Building, Baoan District, Shenzhen, China

Factory Add: 5th Floor, Liutang Industrial Zone, Baoan District, Shenzhen, China

Tel: +86 755 36673336 Fax: +86 755 33130798

Web: www.promed-medical.com Email: info@promed-medical.com

SpO2 Module (Optional)

Measurement Range: 30-100%, ±2% between 80%-90%, others ±5%

Alarm Range: User set high limit and low limit
Alarm Accuracy: ±10% within setting values
Alarm time Accuracy: less than 12 sec

Recorder

Paper: 50mm thermal
Speed: 12.5mm/sec, 25mm/sec, 50mm/sec.
User-selectable 6-second delay
Printing Method: High-resolution, thermal array print head
Print-out Modes: Manual or automatic, user-configurable
On/Off Control: Front panel and paddle
Automatic Function: 9-second recording initiated by alarm activation or defibrillator charge or defibrillator discharge

AED Mode

AED Function: Auto analyze and charge x3 with programmable auto energy level selection, screen prompts, and voice prompts

Shockable Rhythms: Ventricular fibrillation
Ventricular tachycardia
QRS complex wave duration

Charge Control: Control on device front panel, press key on paddle

Manual Mode

Energy Selection: Selectable at 2, 5, 7, 10, 20, 30, 50, 70, 100, 150, 200, 300, 360 joules

Synchronized Mode: Synchronizes defibrillator pulse to patient's R-wave. "SYNC" message displayed on monitor