

DF-300A

Defibrillator



Non-synchronizer : external defibrillator

Sine wave : Monophasic technology

Energy : 0, 20, 50, 100, 160, 250, 300, 360 joule (nominal at a resistance of 50)

Charging time : at 360J less than 10 sec

Paddle options : reusable external adult paddles(pediatric paddles integrated)

Power Supply : operates from the mains(AC)and battery

Battery : charge 3.5h to full ;

full charge to 35 shocks(360 Joules)

sealed rechargeable 12 volt NIMH battery

.Power requirements.

(a) AC Power supply 110V, 60 Hz /220V,50 Hz

(b) Vehicle Voltage DC 12V

Physical characteristics:

Main unit: Portable

Dimensions:L130xW460xH410 Weight:11kg

promed