90111201

BLOOD DONOR CHAIR

DOLSAN

FUNCTIONS



Backrest Adjustment: 0 - 60°



Mattress platform height: 64,5 cm



Trendelenburg position: 0 - 12°

GENERAL SPECIFICATIONS

- Antibacterial and cleanable artificial leather covered, high density foam mattress.
- Backrest position is adjusted between 0-60 degrees by gas spring.
- Height adjustable IV stand.
- Towel rack behind of backrest.
- Trendelenburg position is adjusted between 0-12 degrees by gas spring.
- Armrests height and its angle are adjusted by mechanical system easily.
- Epoxy painted steel construction is resistant corrosion.
- Moves easily by 100 mm diameter four castors two of them with brake to fix the position.



* RV : Vinlex color options page 81.