

## FUNCTIONS



Back plate up : 0 - 75°



Back plate down : 0 - 22°



Trendelenburg position : 0 - 18°



Anti - trendelenburg position : 0 - 18°



Tilt left : 18°



Tilt right : 18°



Leg plate up / down : 10° / 90°  
Head plate up/down: 50° /90°



Leg plate horizontal extension : 90°



Longitudinal slide of tabletop : 15 cm



Longitudinal slide of tabletop : 15 cm





Minimum height : 71 cm



Maximum height : 111 cm

## OPTIONS

### GROUP 1 / ACCESSORY OPTIONS



### DESCRIPTION OF OPTION CODES

**A0** WITHOUT ACCESSORIES

**A1** CONTROL PANEL

**A2** ARM BOARD

**A3** LATERAL SUPPORT

**A4** LEG CRUTCHES

**A5** ANESTHESIA SCREEN FRAME

**A6** RESTRAINT STRAP

Do not select the others if you select without accessory.  
If you select the others you can select all or 1-2-3 ... of them.

## GENERAL SPECIFICATIONS

- In order to attain a high degree of process security in the OR, every movement has to meet a virtually unlimited variety of possibilities for ideal patient positioning in all surgical disciplines, hard and durable stainless steel has to be selected, unique design for two kickout-pin oil cylinder on the base be used and forward-looking technology for stable oil flowing with precise controller undertaken, so that man and machinery could harmonize perfectly.
- Well-designed as well as fine closed steel column ensures hygiene and reduces the cost for cleaning, disinfection and rusty possibility.
- Gas spring activate head and leg plates .
- Special materials for tabletop pave the way for x-ray technology, compatible with interventional radiology.
- Connection with all kinds of accessories ensures different surgeries
- Provide longitudinal slide, guaranteeing effective use of interventional radiology during surgery, especially orthopedic operation.
- Foot-control pedal ensures tabletop's movement when power is off.

## TECHNICAL SPECIFICATIONS

■ Length of tabletop	: 204 cm
■ Width of tabletop	: 52 cm
■ Safe Operating Load (Inc. Mattress and Accessories)	: 180 kg
■ Weight	: approx 150 kg
■ Volume	: 0,68 m <sup>3</sup>
■ Fixation of table	: Mechanically

## ELECTRICAL SPECIFICATIONS

■ Power Supply	: 110~240V 50~60Hz
■ Power	: 850VA

