90103101 EMBRACE 1

BABY CODE

FUNCTIONS



Height: 85,5 cm



Trendelenburg position: 0 - 14°



Anti - trendelenburg position : 0 - 14°

OPTIONS

GROUP 1 / CASTOR OPTIONS





GROUP 2 / CABINET OPTIONS





DESCRIPTION OF OPTION CODES

TO WITHOUT CASTORS

CASTOR TWO WITH BRAKE (50 MM)

D0 WITHOUT CABINET

METAL WITH CABINET

Must select only one of them.

GENERAL SPECIFICATIONS

- Manufacture from antibacterial, transparent flexi glass and portable.
- Trendelenburg and reverse trendelenburg position are adjusted by mechanical system.
- Epoxy painted steel construction is resistant to corrosion.
- Moves easily by 50 mm diameter four castors two of them with brake to fix the position. (Optional)

TECHNICAL SPECIFICATIONS

■ Width : 45 cm
■ Length : 79 cm
■ Weight (Without Accessory) : ~11 kg
■ Weight (Cabinet) : ~8 kg
■ Volume : 0,34 m³





