



Technical and Safety Specification

Crocodile clip

REF 26 600 10

TCE2 - The technical file is created by the manufacturer

Created: 29.03.2017

Update: 29.03.2017

Revision letter: a

Page 1 / 2

- Device classification: Class 1
- Complies with the 93/42/CEE directive

EC REP C.K. Medical Int'l - Laan van Cattenbroeck 10 - 3703BM Zeist - The Netherlands



General description	Crocodile clip	
Attachments	Yes	Photo, diagram
Batch reference	Yes	
ACL Code	/	
Connected device	Yes	Cables and tab electrodes
Single use device	No	
Sterile device	No	
Sterilisation	No	
Autoclavable	No	
Cleaning	Yes	With a dry cloth impregnated with alcohol.
Unsuitable substances		
Component parts	No	
Material		
Labelling	Yes	On the sachet: Denomination, reference, batch, COMED, CE
Packaging	Yes	Per 10 units in a plastic bag
Special packaging	Yes	Temperature between 0 and +50 °C Temperature between 0 and +50 °C
Duration of use	Average	
Marking on medical device	No	
References + extensions	26 600 10	
Designation	Crocodile clip	
Dimensions L x W x H (mm) +- 3mm	57x 20 x 20	
Net weight (g) +- 1%	5	
Colours	White	
Temperature resistance (°C)	0 °C to +50 °C	
Usage	Device for attaching tab electrodes to cables.	
Instruction for use	Push down on the end of the clip to open the "jaw", place it around the electrode tab, then release the pressure so the clip clasps onto the electrode tab	
Precautions for use and/or single use warning	Ensure the device is thoroughly cleaned and disinfected after each use.	
Recycling	No	

Editor: Marion KLEIN
Role: QA
Initials: MK

COMED SAS
8 rue Louise Michel
67200 STRASBOURG - France
www.comed.fr



Technical and Safety Specification

Crocodile clip

REF 26 600 10

TCE2 - The technical file is created by the manufacturer

Created: 29.03.2017

Update: 29.03.2017

Revision letter: a

Page 2 / 2

<i>Recovery</i>	In medical waste containers
<i>Latex/Antistatic</i>	No/No
<i>WEEE/REACH</i>	No/Yes
<i>Phthalates/PVC</i>	No/No
<i>RoHS</i>	NC
<i>Safety data</i>	R-Phrases = None Combination phrases = None

Editor: Marion KLEIN
Role: QA
Initials: MK

COMED SAS
8 rue Louise Michel
67200 STRASBOURG - France
www.comed.fr