



Technical and Safety Specification  
**Mouthpieces for PEAK FLOW METER**  
**REF 30 100 22 - 30 100 30**  
TCE2

Created on: 21.08.2015  
Updated: 21.01.2016  
Index: b  
Page 1/1

- Device classification: Class I
- Complies with the 93/42/CEE directive




Diagramm Halbach - Am Winkelstück - 58239 SCHWERTE - GERMANY

<i>General description</i>	Mouthpieces for Peak Flow Meter	
<i>Attachments</i>	Yes	Photo
<i>Batch Reference</i>	No	
<i>Connected device</i>	Yes	Peak Flow meter
<i>Single use device</i>	Yes	
<i>Sterile device</i>	No	
<i>Sterilisation</i>	No	
<i>Autoclavable</i>	No	
<i>Cleaning</i>	No	
<i>Unsuitable substance</i>	Products that attack cardboard	
<i>Component parts</i>	No	
<i>Material</i>	Cardboard	
<i>Contains latex</i>	No	
<i>RoHS</i>	No	
<i>Labelling</i>	Yes	Name, reference, CE, COMED, address, Batch, Serial no., manufacturer
<i>Packaging</i>	Yes	In a cardboard box
<i>Special packaging</i>	Yes	In a place with an ambient temperature between 0-50°C and relative humidity between 10-95%. In a clean, dry place away from sunlight.
<i>Marking on medical device</i>	No	
<i>References</i>	30 100 22	30 100 30
<i>Designation</i>	Mouthpieces children	Mouthpieces adults
<i>Dimensions L x Ø (mm)</i>	66 x Ø 22	66 x Ø 30
<i>Net weight (g)</i>	3	5
<i>Colours</i>	White	White
<i>Usage</i>	Disposable mouthpieces for PEAK FLOW meter.	
<i>Instruction for use</i>	Place the mouthpiece into the slot on the Peak Flow meter before blowing inside.	
<i>Precautions for use</i>	Disposable mouthpieces. Do not reuse, contamination between patients.	
<i>Recycling</i>	No	
<i>Recycling/Recovery</i>	Waste bin	
<i>WEEE/REACH</i>	No/NC	
<i>PVC/PHthalATES</i>	No/No	
<i>Safety data</i>	R-Phrases = None    Combination phrases = None	

Editor: Marion KLEIN  
Role: QA  
Initials: MK

COMED sas  
8 Rue Louise Michel  
67200 STRASBOURG France  
[www.comed.fr](http://www.comed.fr)