



Safety and Technical Datasheet  
**Measuring tape**  
 REF. 44 447 10 / 44 447 10N  
 TCE3 - Non medical device

Creation : 09.03.2010  
 Mise à jour : 27.10.2016  
 Indice de révision : b

<i>General description</i>	Measuring tape with automatic retraction system by pressure on the button. Length 200 cm or 80 inches.	
<i>Attachments</i>	Yes	Photo
<i>Lot reference</i>	No	
<i>ACL Code</i>	No	
<i>Device connected</i>	No	
<i>Single use device</i>	No	
<i>Sterile device</i>	No	
<i>Sterilization</i>	No	
<i>Autoclavable</i>	No	
<i>Cleaning</i>	Yes	With hand with a product not attacking the ABS
<i>Unsuitable substance(s)</i>	Product attacking the ABS	
<i>Component parts</i>	Yes	Case and tape
<i>Contents</i>	Case made of ABS and tape made of flexible plastic	
<i>Labeling</i>	No	
<i>Packaging</i>	Yes	Individual in plastic bag and case of 20 pieces
<i>Special storage</i>	Yes	At room temperature in a clean and dry place
<i>Duration of use</i>	Temporary	
<i>Marking on the item</i>	COMED Logo	No
<i>References + extensions</i>	44 447 10	44 447 10N
<i>Designation</i>	Measuring tape	
<i>Dimensions L x W x H (mm)+- 3mm</i>	Case = 051 x 13.5 Tape = 2000 x 8 x 0.5	
<i>Net weight (g)+- 1%</i>	29	
<i>Colors</i>	White	
<i>Temperature resistance ( ° C)</i>	-10 to 45° C	
<i>Use</i>	Device used to measure human body parts.	
<i>Instructions for use</i>	Hold the case firmly in one hand and then unroll the tape by gently pulling on the strip that goes beyond. Once the measurement is done, retract the tape in the case by pressing on the button at the center of the case. For proper retraction, ensure that the tape is not twisted.	
<i>Precaution for use and/or single use warning</i>	Clean or disinfect after each contact with a patient.	
<i>Recycling</i>	No	
<i>Disposal</i>	In common waste bins	
<i>Latex/ Antistatic</i>	NC/No	
<i>WEEE/ REACH</i>	No/NC	
<i>Phthalates /PVC</i>	No/No	
<i>RoHs</i>	NC	
<i>Safety data</i>	Phrases R = Null	Combination of phrases = Null

Editor : Marion KLEIN  
 Role : AQ  
 Paraphe : MK

COMED SAS  
 8 rue Louise Michel  
 67200 STRASBOURG - France  
[www.comed.fr](http://www.comed.fr)