

PARC DES FORGES
8 RUE LOUISE MICHEL
67200 STRASBOURG
France
Tel : 03.88.77.13.00
Fax : 03.88.77.00.70
info@comed.fr
www.comed.fr



INSTRUCTIONS

REF 43 653 XX / 43 655 XX / 43 657 XX

SINGLE AND DOUBLE-SIDED STETHOSCOPE

Essential requirement of Directive 93/42/EEC (Annex I Section 13.1):

"The packaging of every device must contain the instructions for use. An exception is allowed for class I and IIa devices, if they can be used safely without the help of such instructions."

For medical examinations. Put the ear tips into your ears and place the diaphragm onto the body part in question. Turn the stethoscope head to listen with the small non-diaphragm side with the non-chill ring (the bell)

PRECAUTIONS FOR USE:

This device should only be used by a competent and qualified person.

Clean the used side with a suitable disinfectant after each examination to avoid contamination.

Ensure the ear tips are clean and free from blockages.

Single-person use, avoid sharing.

Do not fully stretch the tubing as it could break.