

# LS INFRA

LS Infra is an infrared light that relieves the suffering caused by sore muscles and joints. The comfortable warm created by the infrared light penetrates deep into the skin, stimulates blood circulation and warms the muscles. While the muscles are relieved, they automatically stretch and relax. At the same time, the blood circulation helps your body rid itself of impurities and quickly sends oxygen-rich blood to stressed or aching muscles, causing an effective pain relief.

Ease use, height and head orientation can be adjusted. Head 21 centimeters diameter with metal protection grip. Available in 150W or 250W bulb.

It can be adapted to: trolley stand, table clamp (included in standard version), rail support, wall support and extension wall arm.



## LS INFRA

Technical Details	
Light source	IR E-27 150W
	IR E-27 250W
Head diameter	21 cm
Socket	E-27 ceramic
Arm length	90 cm
Colour	White

### LS Infra versions:

- **LS Infra 150** (Cod. E020412)
- **LS Infra 250** (Cod. E02041250)



..... The LS Infra Timer version: 0 - 15 minutes timer integrated on the connection cable.

## LS INFRA TIMER

Technical Details	
Light source	IR E-27 150W
	IR E-27 250W
Head diameter	21 cm
Socket	E-27 ceramic
Arm length	90 cm
Colour	White
Timer 0 -15 minutes	

### LS Infra Timer versions:


- **LS Infra 150 Timer** (Cod. E020412T)
- **LS Infra 250 Timer** (Cod. E02041250T)



**Wall support**  
Cod. B0101800



**AH table clamp**  
Cod. B010190



**Rail Plus bracket**  
Cod. 82026



**Extension wall arm**  
Cod. B0101S0

**Trolley stand**  
4 kg Cod. 09655