

# LUXIFLEX LED

Luxiflex is a cold light luminaire, versatile, of compact design, special for the examination, with great benefits.

## Flexibility

Flexible arm height adjustable 85 cm.  
It offers a wide reaching, maximum comfort and flexibility.

## Compact design

Luxiflex offers easy use, thanks to its compact LED head 8 cm in diameter and the switch located at the base of the lamp, makes this model an ideal luminaire for the most demanding examinations.

## Mounting options


Standard version includes AH table clamp.  
As accessories, we offer a wide selection of supports such as wall clamps, wall extension arm and trolley stand.

## High performance LED

Luxiflex LED 15.000 lux at 50 centimeters.

## LUXIFLEX LED

Technical Details (Cod. I0104AOLED)	
Light source	LED 6W 15°
Lamp lifetime	50.000 hours
Lux a 50 centimeters	15.000 lux
Colour temperature	4.500 °K
Colour rendering	80%
Arm length	72 cm
Colour	White





**Wall support**  
Cod. B0101800



**AH table clamp**  
Cod. B010190



**Rail Plus bracket**  
Cod. 82026



**Extension wall arm**  
Cod. B0101S0



**Trolley stand**  
4 kg Cod. 09655  
8,8 kg Cod. 09651