

# MIMLED

Compact LED Surgery lamp.

Available in two versions: **MIMLED 1000** and **MIMLED 600**.

## HIGH EFFICIENCY

MIMLED 1000                      100.000 lux @ 1 meter

MIMLED 600                        60.000 lux @ 1 meter

## HIGH COLOUR RENDERING

The high colour rendering  $Ra > 94\%$ , describes the exceptional quality of white light, allows tissues to be distinguished with utmost precision.

## LIGHT INTENSITY DIMMING

Adjustable light intensity integrated in the handgrip and on the head which can be regulated using the panel control.

## ENERGY SAVING LIGHT

Low power consumption.

33W      MIMLED 1000

25W      MIMLED 600

## AVAILABLE IN 3 OPTIONS

- Trolley Stand.
- Wall version.
- Ceiling version (up to 3.10 meter height).



MIMLED 1000



MIMLED 600



## MIMLED 1000

Technical Details	
Lux at 1 meter	100.000 lux
Lamp lifetime	50.000 hours
Colour temperature	4.500 °K
Colour rendering	≥ 94%
Light field diameter	195 mm
Power supply	100/230 V - 50/60 Hz
Consumption	33 W

## MIMLED 600


Technical Details	
Lux at 1 meter	60.000 lux
Lamp lifetime	50.000 hours
Colour temperature	4.500 °K
Colour rendering	≥ 94%
Light field diameter	195 mm
Power supply	100/230 V - 50/60 Hz
Consumption	25 W



*Adjustable light intensity  
integrated in the grip*



**Ceiling version**  
Cod. ML1000C  
ML600C



**Wall version**  
Cod. ML1000W  
ML600W



**Trolley Stand version**  
Cod. ML1000FL  
ML600FL

