

VISTA HF

A Woods lamp with magnifier lens. Ultraviolet radiation luminaire specially used in dermatology, cosmetics, veterinary and research.

It provides a polished aluminum reflector. Special magnifying glass made of crystal.

Available versions:

- VISTA HANDLE HF: Hand version with white colour head and black colour grip.
- VISTA PLUS HF: Double arm articulation version.



Protective sliding cover
closed position



VISTA HANDLE HF

Technical Details	(Cod. E051LAT212)
Light source	FL UV compact 9W
Wave length nanometers	360 nanometers
Glass lens	3 diopters
Lens dimensions	10 x 6 cm
Colour	White



Front view and profile view



VISTA PLUS HF

Technical Details	(Cod. E051CAO012)
Light source	FL compact 9W 840
Light source (optional)	FL UV compact 9W
Glass lens	3 diopters
Lens dimensions	10 x 6 cm
Colour	White

VISTA PLUS HF versions:

- **AH table clamp** (Cod. B010190)
- **Wall support** (Cod. B0101800)
- **Rail Plus bracket** (Cod. 82026)
- **Extension wall arm** (Cod. B010150)

