

Afa EX-SERIES STEAM STERILIZER



Afa EX-SERIES

STEAM STERILIZER

ECONOMICAL

SEMI-AUTO

CONTROL

SYSTEM





EX SERIES STEAM STERILIZERS

FEATURES:

1. HORIZONTAL SEMI-AUTOMATIC STERILIZER.
2. WITH DRY AND RE-DRY FUNCTION, AUTO EXHAUST.
3. WITH LOW WATER INDICATOR.
4. WITH AUTO EXHAUST FUNCTION
5. WITH STERILIZATION TEMPERATURE 121C AND 134C BOTH.
6. DOOR CLOSE KNOB IS 90 DEGREE CLOSE/OPEN DESIGN: LIGHT AND CONVENIENT FOR OPERATION.

SPECIFICATIONS:

MODEL	AS-16EX	AS-25EX	AS-40EX
OVERALL	525W x 524D x 425H MM	525W x 691D x 425H MM	575W x 795D x 460H MM
CHAMBER DIAMETER	261MM	261MM	315MM
CHAMBER LENGTH	300MM	470MM	520MM
CHAMBER CAPACITY	16 LITERS	25 LITERS	40 LITERS
STERILIZATION TEMPERATURE/PRESSURE	121 C(1.05KG/CM2) / 134 C (2.07KG/CM2)	121 C(1.05KG/CM2) / 134 C (2.07KG/CM2)	121 C(1.05KG/CM2) / 134 C (2.07KG/CM2)
DESIGN TEMPERATURE	140 C	140 C	140 C
FILLING WATER	MANUAL	MANUAL	MANUAL
DRY FUNCTION	YES	YES	YES
RE-DRY FUNCTION	YES	YES	YES
EXHAUST	AUTO	AUTO	AUTO
POWER CONSUMPTION	1.8KW	2.1KW	3.0KW
VOLTAGE/FREQUENCY	230V/50-60HZ	230V/50-60HZ	230V/50-60HZ
SAFETY DEVICES			
90 DEGREE DOOR OPEN/CLOSE KNOB	YES	YES	NO(SCREW IN TYPE)
OVER TEMPERATURE PROTECTION	YES	YES	YES
OVER PRESSURE PROTECTION	YES	YES	YES
PRESSURE CONTROL SWITCH	YES	YES	YES
ELECTRIC OVERLOAD PROTECTION	NO-FUSE BREAKER	NO-FUSE BREAKER	NO-FUSE BREAKER
CHAMBER WATER LEVEL SENSOR	YES	YES	YES

STANDARD ACCESSORIES:

1. SUS#304 STAINLESS STEEL STERILIZATION BOX x 1
2. SUS#304 STAINLESS STEEL HEATER COVER x 1

OPTIONAL ACCESSORIES:

1. SUS#304 STAINLESS STEEL STERILIZATION BASKET OR
2. TRAY SET (INCLUDING: RACK(SUS304) x1 + TRAY (ALUMINUM) x 3 + PULLING HOLDER(SUS304) x 1)



AFa INDUSTRIES MED LTD.

TEL: 886-2-8512 1168 FAX: 886-2-8512 1116

HTTP: WWW.AFALTD.COM.TW E-MAIL: afa.ltd@msa.hinet.net