

EC-216[®]

universal trolley to
store and transport
flexible endoscopes



Universal
trolley for
trasportation





Vers. 2/2017 gianlucaferrini.com

The **EC-216[®]** equipped with 5 universal fixed shelves stores any brand and models of flexible endoscopes with or without dedicated trays. It is provided with 2 handles to move it easily and 4 medical grade wheels two of which are fixed and 2 swivel ones. Designed with 4 mm thick tempered closing door fixed to the trolley with 2 stainless steel hinges, alternatively with stainless steel closing door and key lock. In case of closing door the leaf has closed-cell silicone gasket.

MODELS

- EC-216[®]**: Trolley without closing door
- EC-216^{®/V}**: Trolley with tempered glass door
- EC-216^{®/I}**: Trolley with stainless steel door and key lock



MEDICAL DEVICES GROUP Srl
 Via G.Fracastoro, 1 48124 Ravenna (RA)
 Tel +39 0544 1826013 - Fax +39 0544 1826015
 www.mdg-srl.com - info@mdg-srl.com



Directive 93/42/EEC Annex VII compliant