



Handheld Urine Analyzer



Convenient

Small in size and light in weight



Easy

Onscreen display or easy printing via a handheld printer



Efficient

Automatically captures and reports strip results



Fast

Only 70 seconds per strip



Reliable

High reproducibility, semi-quantitative urine chemistry



Application

Community Hospitals, Hospitals, Public Health Centers, Health Care Centers, Home Care

11 Test Items

- Leukocytes
- Urobilinogen
- Nitrite
- Protein
- PH
- Blood
- Specific Gravity
- Ketone
- Bilirubin
- Glucose
- Vitamin C

12 Test Items: 11 Test Items + Microalbumin




14 Test Items: 11 Test Items + Microalbumin Creatinine Calcium

Wireless printer

- optional



Specification Comparison

	Ui Series 	Ui-2A Series 	Ui-10A Series 
Screen	LCD screen	3.5" TFT+Touch Screen	No screen
Key pad	Capacitive touch key		
Speed	140 tests / hour (fast mode), 50 tests / hour (normal mode)		
Test items	11	11/12/14	
Dimension	110*68*27mm	106*63*27.5mm	110*62*27.5mm
Capacity	1000 recent test results		
Printer	Wireless Thermal Printer (optional)		
Interface	Mini USB	Micro USB	
Battery	Lithium battery		AAA Dry battery
Bluetooth	✓		
WiFi	✓		

CHENGDU EMPSUN MEDICAL TECHNOLOGY CO., LTD.

Add: No.333, Hezuo Road, West Hi-tech Zone, Chengdu, Sichuan, China.

Domestic Sales Tel: +86-28-6758-5399

International Sales Tel: +86-28-6558-9558

Customers Service Tel: +86-400-605-9899

E-mail: export@empsun.com

Website: www.empsun.com/en/index.php

All rights are reserved by Chengdu Empsun Medical Technology Co., LTD. like changing product item number, specification, appearance, etc. Pls be subject to the latest leaflet. EN.Ui.201902



Official WeChat