



HF SERIES

BREATH RESUSCITATION BAG

Stock code [300314]



Brief Introduction

Main parts are made of silicon and PC material;
Aspirator valve is equipped with the fish mouth shape valve which can be shut automatically while patient is breathing out and prevent reverse flow;
Equipped with the pressure-limiting valve to avoid the higher pressure output.

Standard Configuration

Air Bag, Aspirate Valve, Air Collection Bag, Mask, Fish Mouth Shape Valve.

Specification

Scope: HF-I is for the patient's weight < 5kg
HF-II is for the patient's weight < 10kg
HF-III is for the patient's weight > 40kg
Output pressure: 3~4.5Kpa (HF-I/HF-II)
≤ 6Kpa(HF-III)
Tidal volume: ≥ 50ml(HF-I)
≥ 150ml(HF-II)
≥ 600ml(HF-III)
Ventilation frequency: 12 times/min, manually control
(HF-I/HF-II/HF-III)

Product Packaging

Package: The whole unit is packed into one carton
Dimension: 260mm x 200mm x 110mm (HF-I/HF-II)
350mm x 250mm x 170mm (HF-III)
Gross weight: 0.8kg (HF-I / HF-II)
1.0kg (HF-III)



HF-I



HF-II



HF-III



Ningbo David Medical Device Co.,Ltd

Add: Bldg No. 5 Ningbo Smart Park, #98 Chuangyuan Rd, Hi-tech Zone, Ningbo, China
Tel: 0086-574-8780 0008 or 8780 0007 /Fax: 0086-574-8780 3714
E-mail: sales@nbdavid.com /Web.: www.nbdavid.com /P.C.: 315048