

90116101

IV STAND

GENERAL SPECIFICATIONS

- Height of IV stand is adjusted owing to mechanic system
- Long lived metallic body and three legs bottom epoxy painted steel construction
- There are two serum hooks.
- Moves easily with 50 mm diameter castor.

TECHNICAL SPECIFICATIONS

- Height (min) : 107 cm
- Height (max) : 192 cm
- Weight : approx 2,5 kg
- Volume : 0,23 m³



90116201

IV STAND

GENERAL SPECIFICATIONS

- Height of IV stand is adjusted owing to mechanic system.
- Long lived plastics five legs bottom.
- There are four serum hooks.
- Epoxy painted steel construction is resistant to corrosion.
- Moves easily with 50 mm diameter castor.

TECHNICAL SPECIFICATIONS

- Height (min) : 110 cm
- Height (max) : 195 cm
- Weight : approx 3,5 kg
- Volume : 0,33 m³

