



- 11 parameters
- 120 tests/hr
- The color and limpidness can be input
- Original matching urine strips for best accuracy

KH-100

Urine Analyzer

Jinan Kinghawk Technology Co.,Ltd

KH-100 Urine Analyzer

Performance

- LCD monitor
- Adopt luminescence source with high brightness to minimize interference of the environmental light.
- Date and patient number can be input.
- Single test and continuous test are available
- All the test data can be saved automatically.

Technical Specifications

- **Wave length:** 440- 660nm
- **Throughput:** 120 test/hr
- **Storage:** 1000 patient data
- **Display:** 240x64 LCD monitor
- **Printer:** built- in thermal printer
- **Interface:** RS232 serial
- **Calibration:** fully automatic calibration
- **Condition:** temperature: 0- 40°C , humidity:5- 90%
- **Power supply:** 220VAC±10% 50~60Hz 30W
- **Dimension:** 390(L)x320(W)x140(H)mm
- **Weight:** 5Kg(net), 6.5Kg(gross)

Test Parameters

- PH, Nitrite, Protein
- Glucose, occult blood, bilirubin
- Urobilinogen, ketone
- Specific gravity, Leukocytes
- Ascorbic Acid
- Color
- Limpidness

Urine Strips



Jinan Kinghawk Technology Co.,Ltd

No.205, E South Building,
Environmental Protection Technology Park,
Jinan 250101, Shandong, China
Tel: 0086 531 8970 2821/2/3
Fax: 0086 531 8268 7770
E-mail: export@kinghawktech.com

Profession & Care