

## The First Generation - LG101 (G1)



70 minutes battery life, 5 liters continuous flow,  
**Non-Medical grade**

### Economical Portable Oxygen Concentrator

LG101 could use both in car and by battery. It produce 40% to 93% oxygen with 5 liters continuous flow. It offer customers another option which has limited budget

2 0 1 3

## Features

- 93% @ 1LPM
- 40% @ 5LPM for healthcare use
- China first POC with continuous flow
- Use at home / in the car / for travel



### Suitable for:

Mild disease patients  
Normal people health care  
Plateau tourism

Model: LG101 (G1)

Weight: 12 lbs (5.4 kg)

Size: 10.2in x 7.7 in x 15.2 in  
26 cm x 19 cm x 38 cm


Warranty	2 Years for machine
	1 Year for sieves

### Flow settings and oxygen saturation

Continuous flow settings					
Continuous Flow	1LPM	2LPM	3LPM	4LPM	5LPM
Oxygen Saturation	93%	75%	55%	45%	40%
Battery Life	70 Mins	70 Mins	70 Mins	70 Mins	70 Mins

### Working Environment

Working Altitude	up to 12,000 ft (3657 meters) support use in plateau, mountains and high altitudes
Working Sound Level	≤ 40 dB (nominal) at flat table, ≤45dB when travelling
Operational Temperature	41° F to 104° F (5° C to 40° C)
Relative Humidity	15% to 95%



Demonstration Video



Package information	
Inner color box	54*23*30 cm
Outside box (pack 2 pieces)	55*48*30 cm
20 GP container	560 pieces



Accessories	Photo	Item Code	Quantity
Battery		1608	1 piece
Charger		1609	1 piece
Car Inverter		1610	1 piece
Spare Air Filter		1703	1 piece
Nasal Cannula(2 meter)		1704	1 piece
Humidifier Bottle		1705	1 piece
Shoulder Bag		1902	1 piece
Trolley Cart		1903	1 piece
User Manual		1904	1 piece

