

LUPA LED HF

LED magnifier with round lens, HF. LED Magnifier is intended to perform a meticulous work and with great accuracy. Its use is ideal for hospitals, measured studies, laboratories, aesthetic clinics, etc.

Functionality and Performance

Generates a light without shadows.
Internal spring balanced arms that allow different positions.
Switch on the head.
Round head with convex lens, adjustable and simple mechanism

Mounting Options

Standard version includes AH table clamp.

As accessories, we offer a wide selection of supports such as wall clamps, wall extension arm and trolley stand.

Cleaning and Safety

Head exempt of heat.
Protective cover on the head with hinge mechanism to prevent from dust and scratches that may damage the lens.

Economy

High performance, low power 6W.
Long life, 40,000 hours.
Durable materials.
Very low maintenance costs.



LUPA LED HF

Technical Details		(Cod. 09340)
Light source	LED HF 6W	
Lamp lifetime	40.000 horas	
Colour temperature	6.500 °K	
Glass lens	5 diopters	
Lens diameter	127 mm	
Arm length	90 cm	
Weight	3 kg	
Colour: White	White	



Trolley stand
8,8 kg Cod. 09651






Wall support
Cod. B0101800




AH table clamp
Cod. B010190



Rail Plus bracket
Cod. 82026




Extension wall arm
Cod. B010150




Fixing lock

Hinge mechanism



Safety lid

Switch on the head



60 LEDs

Circular glass lens