

LUXIFLEX

Luxiflex is a versatile and stylish lamp for examination and minor surgery procedures, ideal in health centers, clinics and hospitals. It is also suitable for veterinarian examinations and procedures.

The 85 cm arm's length allows an ample turn around providing maximum comfort and flexibility.

Low temperature lamp without infrared radiation. Switch located on the arm's base.

It can be adapted to: trolley stand, table clamp (included in standard version), rail support, wall support and extension wall arm.

Also available in LED version: Luxiflex LED and Luxiflex LED Plus.



LUXIFLEX HALOGEN

Technical Details		(Cod. I0104AO)
Light source	Hal MasterLine ES 35W 8°	
Lamp lifetime	5.000 hours	
Lux a 50 centimeters	50.000 lux	
Colour temperature	4.000 °K	
Colour rendering	95%	
Arm length	72 cm	
Colour	White	
Glass filter for colour temperature		



Wall support
Cod. B0101800



AH table clamp
Cod. B010190



Rail Plus bracket
Cod. 82026



Extension wall arm
Cod. B0101S0



Trolley stand
4 kg Cod. 09655
8,8 kg Cod. 09651