

## 90104101

## MATTRESS

**GENERAL SPECIFICATIONS**

- The foam is for hospital beds use.
- The cover is non latex and washable.
- The cover is conductive, antistatic, antibacterial, antiacid,antifungus.
- The cover is waterproof and easy to clean.
- The cover is non latex and washable.

**TECHNICAL SPECIFICATIONS**

■ Width	: 90 cm
■ Depth	: 200 cm
■ Height	: 15 cm
■ Volume	: 0,24 m <sup>3</sup>



## 90104103

## MATTRESS VISCO

**GENERAL SPECIFICATIONS**

- Special foam sponge is able to be take a shape according to the hospital bed positions.
- Anti-bacterial bed cover provides air circulation and protects the bed with its liquid tight feature.
- Bed cover, which can be removable and assembled by a zip fastener system if necessary, is of washable feature.

**TECHNICAL SPECIFICATIONS**

■ Width	: 90 cm
■ Depth	: 200 cm
■ Height	: 15 cm
■ Weight	: ~ 8 kg
■ Volume	: 0,21 m <sup>3</sup>

