



Microplate Washer W206



Sweep Function
Aspiration



1~99 adjustable washing
cycles



Automatic Liquid
Monitoring System



Multi-buffer Module



8/12 type washing
head-compatible with F,U or
V-bottom Plates



Hospital, Quarantine, Lab,
Clinic, Veterinary

Details



Microplate

Technical Specifications

Plate type	96-well plate	Temperature	5°C~35°C
Washing head	8/12 type	Humidity	15%~85%
Residual aspiration volume	< 2ul		
Dispensing volume	10~3000ul, step by 1ul		
No. of washing cycle	1~99 adjustable		
Soak time	0~999s		
Shaking time	0~999s		
Weight	8Kg		
Dimension (mm)	400(L) x 310(W) x 169(H)		
Power supply	AC 100V~240V, 50/60Hz		

CHENGDU EMFSUN MEDICAL TECHNOLOGY CO., LTD.

Add: No.333,Hezuo Road,West Hi-tech Zone,Chengdu,Sichuan,China.

Domestic Sales Tel: +86-28-6758-5399

International Sales Tel: +86-28-6558-9558

Customers Service Tel: +86-400-605-9899

E-mail: export@empusun.com

Website: www.empusun.com/en/index.php

All rights are reserved by Chengdu Empsun Medical Technology Co.,LTD.like changing product item number,specification,appearance,etc. Pls be subject to the latest leaflet.EN.W206.201902



Official WeChat