## 90110201 NEGATOSCOPE

# 90110202 NEGATOSCOPE

# DOLSAN

#### **GENERAL SPECIFICATIONS**

- Screen is made of white plexiglass.
- One film is able to be attached on the screen.
- Screen is designed for removing or attaching film easily.
- 32 Watt round fluorescent is used in the system.
- Steel construction is epoxy painted, which is durable against corrosion.

#### **TECHNICAL SPECIFICATIONS**

: 40,5 cm ■ Width ■ Depth : 10 cm ■ Height : 49,5 cm ■ Weight : approx 6 kg ■ Volume : 0,03 m<sup>3</sup>



#### **GENERAL SPECIFICATIONS**

- Screen is made of white plexiglass.
- Two films are able to be attached on the screen.
- Screen is designed for removing or attaching film easily.
- 32 Watt round fluorescent is used in the system.
- Steel construction is epoxy painted, which is durable against corrosion.

#### TECHNICAL SPECIFICATIONS

■ Width : 75 cm  ■ Depth : 10 cm  ■ Height : 49,5 cm  ■ Weight : approx 6,5 kg  ■ Volume : 0,05 m³	I E O I II II O I E O I E O	011 107 1110110
■ Height : 49,5 cm ■ Weight : approx 6,5 kg	■ Width	: 75 cm
■ Weight : approx 6,5 kg	■ Depth	: 10 cm
11 . 0	■ Height	: 49,5 cm
■ Volume : 0,05 m³	■ Weight	: approx 6,5 kg
	■ Volume	: 0,05 m <sup>3</sup>



## 90110203 NEGATOSCOPE

#### **GENERAL SPECIFICATIONS**

- Screen is made of white plexiglass.
- Three films are able to be attached on the screen.
- Screen is designed for removing or attaching film easily.
- 32 Watt round fluorescent is used in the system.
- Steel construction is epoxy painted, which is durable against corrosion.

#### TECHNICAL SPECIFICATIONS

Width	:	110 cm
Depth	:	10 cm
Height	:	49,5 cm
Weight	:	approx 7,5 kg
Volume	:	0,07 m <sup>3</sup>

