

90106101

OVER BED TABLE MECHANICAL

GENERAL SPECIFICATIONS

- Laminate or ABS plastic top is durable against to water and heat.
- Height of over bed table is adjusted by mechanical handle according to height of bed.
- Epoxy painted steel construction is resistant to corrosion.
- With four castors (two of them with brake) can move freely.

TECHNICAL SPECIFICATIONS





■ Width	: 80 cm
■ Depth	: 40 cm
■ Height (min)	: 79 cm
■ Height (max)	: 103 cm
■ Weight	: ~ 9 kg
■ Volume	: 0,30 m ³

* RS : Melamine color options page 81.

* RL : Laminate color options page 80.

OPTIONS

GROUP 1 / TABLE TOP OPTIONS



DESCRIPTION OF OPTION CODES

- YL1 LAMINATE COATED CHIPBOARD TOP
 - YL2 ABS TOP
 - YL3 COMPACT LAMINATE TOP
 - YL4 MELAMINE COATED CHIPBOARD TOP
 - YL5 TILTED COMPACT LAMINATE UPPER TABLE
- Must select only one of them.



90106201

OVER BED TABLE WITH GAS SPRING

GENERAL SPECIFICATIONS

- Laminate or ABS plastic top is durable against to water and heat.
- High of over bed table is adjusted by gas spring according to height of bed.
- Epoxy painted steel construction is resistant to corrosion.
- With four castors (two of them with brake) can move freely.

TECHNICAL SPECIFICATIONS





■ Width	: 80 cm
■ Depth	: 40 cm
■ Height (min)	: 79 cm
■ Height (max)	: 103 cm
■ Weight	: ~ 10,5 kg
■ Volume	: 0,30 m ³

* RS : Melamine color options page 81.

* RL : Laminate color options page 80.

OPTIONS

GROUP 1 / TABLE TOP OPTIONS



DESCRIPTION OF OPTION CODES

- YL1 LAMINATE COATED CHIPBOARD TOP
 - YL2 ABS TOP
 - YL3 COMPACT LAMINATE TOP
 - YL4 MELAMINE COATED CHIPBOARD TOP
 - YL5 TILTED COMPACT LAMINATE UPPER TABLE
- Must select only one of them.

