

Plug-in Type Patient Monitor

MODULAR PATIENT MONITOR

promed



PM-8 Plus

PM-10 Plus

PM-12 Plus

PM-15 Plus

Features

- Full touch screen selectable, more convenient for operation.
- Real-time S-T segment analysis, pace-maker detection.
- 3-level audio/ visual alarms.
- Efficient resistance to interference of defibrillator and HF knife.
- Lasting working capability with built-in rechargeable lithium battery.
- 400 groups of NIBP list, 6000 seconds ECG waveform recall, 60 alarm event records recall, 7-day trend chart in storage.
- Wire/ wireless central monitoring system applicable to ICU/ CCU/ OR.
- Support solution: wall mounting, mobile trolley.

Standard Configuration

SpO2, PR, NIBP, ECG, RESP, TEMP

Optional Configuration:

Nellcor SpO2, 2-IBP, ETCO2, Touch screen, Thermal printer

Sales Add: Room 705, Jinfulai Building, Baoan District, Shenzhen, China

Factory Add: 5th Floor, Liutang Industrial Zone, Baoan District, Shenzhen, China

Tel: +86 755 36673336 Fax: +86 755 33130798

Http://www.promed-medical.com E-mail: info@promed-medical.com

Introduction

Characterized in half plug-in structure, modular patient monitor dedicates to provide comprehensive health care solutions with more selections for customers. With its plug-and-play modules, modular patient monitor is able to inspect a wide range of parameters simultaneously, adapting from general monitoring to critical care. Additionally, monitoring capabilities such as ETCO2 and 2-IBP are optional.

