

NOVIMED 4 Drawers Hospital Cart (1-0-3)

Reference: F-DNOV-1-0-3



Hospital Carts New Series NOVIMED

New design with innovative advantages and new characteristics designed to make the work of the healthcare staff as easy as possible.

- Lightweighth and maoeuvrable
- More resistant
- Versatil
- Easier to clean
- Best quality price ratio

New Design with:

- Worktop with three raised sides (ABS)
- Sliding table top extension
- 4 Drawers (ABS) with front-label
- DIN side rail
- Bumpers of plastic
- Phenolic base
- 4 Manner wheels of 125 mm, 2 with locking
- One push handle as standar (ABS)
- Anodised aluminum profiles to assembly accessories
- Ball bearin rails total extraction and soft-close system
- New system DIN side rail, easy assemble on aluminium profiles
- Accesories not included

