

## NOVIMED 6 Drawers Hospital Cart (2-4-0)

Reference: F-DNOV-2-4-0



### Hospital Carts New Series NOVIMED

**New design with innovative advantages and new characteristics designed to make the work of the healthcare staff as easy as possible.**

- Lightweighth and maoeuvrable
- More resistant
- Versatil
- Easier to clean
- Best quality price ratio

-----  
New Design with:

- Worktop with three raised sides (ABS)
- Sliding table top extension
- 6 Drawers (ABS) with front-label
- DIN side rail
- Bumpers of plastic
- Phenolic base
- 4 Manner wheels of 125 mm, 2 with locking
- One push handle as standar (ABS)
- Anodised aluminum profiles to assembly accessories
- Ball bearin rails total extraction and soft-close system
- New system DIN side rail, easy assemble on aluminium profiles
- Accesories not included

