90111101 **BLOOD RECEPTION CHAIR**

GENERAL SPECIFICATIONS

- Antibacterial and cleanable artificial leather covered high density foam cover.
- Height and back part chair adjustable by gas spring.
- Height and angle of arm holder adjustable by mechanical system.
- Epoxy painted steel construction is resistant to corrosion.
- Moves easily by 50 mm diameter castors.
- Vinlex colors are optional.(RVX)

TECHNICAL SPECIFICATIONS



90111102 **BLOOD RECEPTION CHAIR**

GENERAL SPECIFICATIONS

- Antibacterial and cleanable artificial leather covered high density foam cover.
- Backwards at the same time moving the foot platform.
- Height and angle of arm holder adjustable by mechanical system.
- Height adjustable IV stand.
- Epoxy painted steel construction is resistant to corrosion.
- Optional artifical leather choice. (RDX)

Width Depth	: 70,5 cm : 96 - 170 cm
Sitting Height	
Backrest Adjustment	: 44 cm : 25°-65°
Weight Volume	: approx 33 kg : 0,90 m ³









