

Company Profile

Since its establishment, SonoScape has been committed to providing high quality medical equipment for the healthcare sector. SonoScape specializes in the development and production of diagnostic ultrasound solutions. By introducing advanced imaging techniques, SonoScape has improved diagnostic accuracy and therefore enabled better health outcomes. Since 2002, SonoScape has reached and benefited millions of people. With the world in mind, SonoScape will continue providing more effective and accessible healthcare solutions, through persistent innovation and passion for life.

Quote from the Analyst of

"The Frost & Sullivan Award for Growth Leadership is presented to the company that has demonstrated excellence in capturing the highest annual compound growth rate for the last 3 years...The sales of ultrasound devices contribute an overwhelming proportion of total revenue, helping the company strengthen its leading position to outperform the competitors in the marketplace. With increasing penetration of mid-to high-end ultrasound devices, SonoScape is expected to see significant growth potential in the future..."

Company Milestone

- 2002: Company Founded in Shenzhen, China
- 2003: Released SSI-1000: the 1st 15" Portable Color Doppler system in China
- 2004: Released SSI-2000: the 1st PC platform Color Doppler system in China
- 2005: Received the "High Technology Company" award from the PRC government
- 2007: Received "CHINA TOP BRAND" award in the Medical Equipment Industry
- 2007: Released the 1st Real time 4D ultrasound system in China
- 2008: Received "European Entrepreneurial Company 2008" award from FROST & SULLIVAN
- 2008: Received "Flagship Company" award in the Medical Equipment Industry in China
- 2009: Received "Product Quality Leadership Award 2009" from FROST & SULLIVAN
- 2011: Received the Reddot 2011 Product Design Award for S20 in Essen, Germany
- 2013: Received "Ultrasound Market Growth Leadership Award, 2013" from FROST & SULLIVAN
- 2014: Received the iF product design award 2014 for S9 in Munich, Germany
- 2014: Received "Company of the Year in Ultrasound Market, 2014" from FROST & SULLIVAN



ISO 13485

CE 0197

SonoScape

Yizhe Building, Yuquan Road, Shenzhen, 518051, China
Tel: 86-755-26722890 Fax: 86-755-26722850
E-mail:sonoscape@sonoscape.net www.sonoscape.com

SonoScape

A Value Choice
Beyond Your Expectation



S11

U-S1120150120

Caring for Life through Innovation

The image shows a Philips S11 ultrasound machine. The main screen displays a color Doppler ultrasound image of a vessel. The control panel features a keyboard, a trackball, and several buttons with green labels. The machine is white and black, with a large monitor on a stand. The background is a light blue and white gradient with a circular pattern of thin lines.

S11

*A Value Choice
Beyond Your Expectation*

Integrated with SonoScape's new software platform, together with its ergonomic design and optimized workflow, S11 will help you diagnose quickly and easily.

SonoScape's speckle reduction technology will enhance the image acquisition and increase diagnostic accuracy. The new patient file management software will enhance the operational workflow. The user-definable pre-sets and shortcut keys will help you work in your own style. The panel and keyboard layout are designed to allow greater functionality for the operator. The powerful new platform will provide a seamless experience for all users.

S11 encompasses value, features and ergonomics.

Various Diagnostic Applications

- General Applications: Abdomen, Obstetrics, Gynecology, Small parts, Muscle-skeletal, rheumatology, Vascular, Urology, Cardiology
- Support Convex, Micro-convex, Transvaginal, Transrectal, Linear, Intraoperative, Phased array, 4D Volumetric transducers

Reliable Color Doppler System

- Imaging modes: B, M, Color, Power, PW, CW(optional)
- Main technologies: Spatial Compound Imaging, Panoramic Imaging, Anatomic M Mode, Freehand 3D Imaging, 4D Imaging (some of them are optional.)



Professional OB/GYN Clinical Solutions

Complete Scanning modes, outstanding 2D performance, sensitive blood flow, freehand 3D and volume 4D imaging technology, endow S11 with the best package for all applications, especially women's healthcare in its class.

Excellent 2D image with advanced transducer techniques

- Quad B mode, convenient AFI measurement
- Endocavity transducer with wide scanning angle and temperature detection
- Specially designed biopsy guide for wide range of transducers
- 3D,4D acquisition enable strong emotion connection for an early experience of motherhood

Specialized OB/GYN software package enables better quantification and calculation

- Obstetrical: GS, CRL, BPD, HC, AC, FL, CER, OFD, Fibula, AA, etc.
- Gynecological: UT-L, UT-H, UT-W, Cx, EN-T, AFI, OV L/H/W, etc.
- Measurement Methods: Hadlock, Jeanty, Kurtz, Hansmann, Nelson, Tokyo, etc.
- EFW Measurement Methods: WEI/SAB, Shepard, Hadlock, Hansmann, etc.





- Smart patient file management: Easy to retrieve patient image and data at any time;
- M-Tuning: one button image optimization for 2D and spectral Doppler
- Complete DICOM system: Save, Storage Commitment, MPPS, Print, Worklist, advanced file transmission/management;
- Versatile connectivity: 320G Hard drive, USB ports, VGA, S-video, DVD RW, LAN port, thermal printer and external USB laser/jet printer support;
- Flexible user-defined functions, integrated function keys, icons, measurement, report, comment, layout...

Exceptional Workflow

S11 offers comprehensive data storage, transmission solutions and intelligent data management. It also offers you easy access to customize the settings according to your own preference.

WORKFLOW OWN





- Smart patient file management: Easy to retrieve patient image and data at any time;
- M-Tuning: one button image optimization for 2D and spectral Doppler
- Complete DICOM system: Save, Storage Commitment, MPPS, Print, Worklist, advanced file transmission/management;
- Versatile connectivity: 320G Hard drive, USB ports, VGA, S-video, DVD RW, LAN port, thermal printer and external USB laser/jet printer support;
- Flexible user-defined functions, integrated function keys, icons, measurement, report, comment, layout...

Exceptional Workflow

S11 offers comprehensive data storage, transmission solutions and intelligent data management. It also offers you easy access to customize the settings according to your own preference.

