

## 90115101

## SCREEN

**GENERAL SPECIFICATIONS**

- Two winged screen is adjusted by casters.
- The curtains are removable, attachable and washable.
- Epoxy painted steel construction is resistant to corrosion.
- Screen is able to be used in narrow area, owing to foldable wings system if necessary.

**TECHNICAL SPECIFICATIONS**

- |          |                       |
|----------|-----------------------|
| ■ Length | : 140 cm              |
| ■ Height | : 166 cm              |
| ■ Weight | : approx 7,4 kg       |
| ■ Volume | : 0,09 m <sup>3</sup> |



## 90115102

## SCREEN

**GENERAL SPECIFICATIONS**

- Three winged screen is adjusted by castors.
- The curtain is removable, attachable and washable.
- Epoxy painted steel construction is resistant to corrosion.
- Screen is able to be used in narrow area, owing to foldable wings system if necessary.

**TECHNICAL SPECIFICATIONS**

- |          |                       |
|----------|-----------------------|
| ■ Length | : 210 cm              |
| ■ Height | : 166 cm              |
| ■ Weight | : approx 10,8 kg      |
| ■ Volume | : 0,12 m <sup>3</sup> |

