## 90101134 SIMPLE 1

### **FUNCTIONS**



Backrest adjustment: 0 - 70°



Thigh rest adjustment: 0 - 40°



Foot rest mechanical adjustment : 0 - 25°





Fowler position



Trendelenburg position: 0 - 12°

## **HOSPITAL BED 2 MOTORS**

### **OPTIONS**



#### **GROUP 3 / OTHER ACCESSORY OPTIONS**









WITHOUT OTHER ACCESSORIES

LIFTING POLE : BEARING CAPACITY 75 Kg

TELEVISION STAND FOR 17" LCD TV

TRACTION APPARATUS



### **DESCRIPTION OF OPTION CODES**

- TWIN CASTOR TWO WITH BRAKE (120 MM)
- CASTOR TWO WITH BRAKE (125 MM)
- METAL CASTOR TWO WITH BRAKE (125 MM)
- WITH METAL SIDE RAIL
- WITH METAL SIDE RAIL
- Must select only one of them.
- Do not select the others if you select without accessory.

  If you select the others you can select all or 1-2-3 ... of them.

V STAND : ONE HEIGHT ADJUSTMENT SYSTEM

### **GENERAL SPECIFICATIONS**

- Easily cleanable metal bed platform designed for the ventilation of the bed.
- Head and foot movements with the help of an electrical motor.
- Trendelenburg movement by the help of shock absorber mechanism at foot side.
- Ergonomic, durable and easily cleanable metal rails system with maximum safety.
- Melamine coated chipboard head and foot panels.
- Removable and installable head system that can be also used as a resultation board if necessary.
- Shock resistance plastic bumper wheels.
- Castor system with two-castor brake in a diameter of 125 mm.
- Metal components dyed with anti-corrosion electrostatic powder paint after being treated with phosphate.
- Optional melamine coated chipboard color choice. (RSX)

#### **TECHNICAL SPECIFICATIONS**

■ External dimension	:	97 x 207 x 50 cm
■ Recommended mattress dimension	:	85 x 190 x 12 cm
■ Safe Operating Load		
(Inc. Mattress and Accessories)	:	230 kg
(Inc. Mattress and Accessories)  ■ Weight		230 kg ~ 86 kg

#### **ELECTRICAL SPECIFICATIONS**

■ Power Input	:	230 V AC 50 / 60 Hz 50W
■ Power Output	:	24 V DC







# **DOLSAN**

