90114304

EXAMINATION COUCH

GENERAL SPECIFICATIONS

- Antibacterial and cleanable artificial leather covered, high density foam mattress with its pillow.
- Backrest adjustment by mechanical system between 0-48 degrees.
 Owing to foldable foot, examination couch is able to be preserved
- and carried with very little volumes.
- Epoxy painted steel construction is resistant to corrosion.
- Vinlex colors are optional.^(RVX)

TECHNICAL SPECIFICATIONS

 Width 	:	59 cm
Length	:	189 cm
 Height 	:	66 cm
 Backrest Adjustment (max) 	:	48°
 Weight 	:	approx 29 kg
Volume	:	0,28 m ³



90118101

STEP

GENERAL SPECIFICATIONS

Easy cleanable laminated stair is resistant to water and heat.

- In order to protect ground against to tearing plastic cover used on the legs.
- Epoxy painted steel construction is resistant to corrosion.
- Optional laminate color choice.^(RLX)

TECHNICAL SPECIFICATIONS

 Width 	: 44 cm	
 Depth 	: 47,5 cm	
 Height 	: 50 cm	
 Weight 	: approx 6 k	g
 Volume 	: 0,11 m ³	



* RL : Laminate color options page 80.

