

90107101

STRETCHER

FUNCTIONS



Backrest Adjustment : 0 - 45°

OPTIONS

GROUP 1 / IV STAND OPTIONS



DESCRIPTION OF OPTION CODES

W0 WITHOUT OTHER ACCESSORIES

W8 IV STAND : ONE HEIGHT ADJUSTMENT SYSTEM

Must select only one of them.

GENERAL SPECIFICATIONS

- Antibacterial and cleanable artificial leather covered, high density foam mattress.
- Backrest position is adjusted between 0-45 degrees by mechanical system.
- Ergonomic, durable and easily cleanable metal rails system with maximum safety.
- Height adjustable IV stand is able to be located on both head and foot side corners.
- There are four bumpers on corners of bed for protection.
- Epoxy painted steel construction is, resistant to corrosion.
- Upper platform is portable in case of emergency, it can be easily removed and transferred from trolley.
- Moves easily by 150 mm diameter four castors, two of them with brake to fix the position.
- Optional vinlex color choice. (RVX)

TECHNICAL SPECIFICATIONS

■ Width	: 71 cm
■ Length	: 200,5 cm
■ Height	: 78 cm
■ Weight	: approx 50 kg
■ Volume	: 1,20 m ³



* RV : Vinlex color options page 81.