

# Temp Sitter Wireless Thermometer



## Parameter

Dimension	51mm*32mm*7mm
Waterproof Class	IPX7
Battery Life	365 days (8 hour/day)
Blue tooth version	4.0
Distance (Blue tooth)	5m
Measurement range	25.00°C~45.00°C
Accuracy	±0.05°C ( When the temp is between 35.00°C~38.50°C) ±0.1°C ( When the temp is between 25.00°C~34.99°C)
Data storage time	10days
Automatically link (blue tooth)	Support
APP platform	IOS, Android
Power on/off	Automatically
Alarm	Sound and message
Cloud storage of data	Support
Family account management	Log in and manage family account



# Temp Sitter

## Wireless Thermometer



**BIOLIGHT CO.,LTD.**

Tel: +86-756-3399935 Fax: +86-756-3399911 E-mail: Overseas@blt.com.cn  
Address: No.2 Innovation First Road, Technical Innovation Coast, Hi-tech Zone, Zhuhai, P.R.China Postcode:519085



# Temp Sitter Wireless Thermometer

# Temp Sitter Wireless Thermometer



Baby getting a fever is always tricky for parents. Let Temp Sitter set you free from the tiresome work of frequent measurement and the worries of hyperthermia harm.

## Traditional Thermometer VS. Temp Sitter



Resistance caused by uncomfortable feelings

Well-behaved owe to the cosy feeling




Unable to monitor continuously

24h continuous monitoring



No in time medical reminder

Automatic high-temp alarm



Toxic and unreliable materials

Materials of safe and high quality



Turn on the APP and start the intelligent temperature measurement.



Cloud Management



Low Power Consumption



High Precision



Water Proof



High Temp Alarm



7 Days Trend Review



Data storage



Long Distance Monitor