# 90108101 TROLLEY DRESSING

## 90108201 TROLLEY DRUG DISTRIBUTION



### **GENERAL SPECIFICATIONS**

- Easy cleanable and resistant ABS plastics top platform.
- At both sides there are one piece of quartet drug boxes produced from transparent plexiglass.
- Sliding shelf and drawer for drug preparation.
- There is serum containers made of ABS plastics.
- There is a stainless steel dustbin and dressing bowl.
- Height adjustable IV stand.
- Epoxy painted steel construction is, resistant to corrosion.
- Moves easily by 100 mm diameter four castors two of them with brake to fix the position.

#### **TECHNICAL SPECIFICATIONS**

| <ul> <li>Width</li> </ul>  | : 93 cm  | ı       |
|----------------------------|----------|---------|
| <ul> <li>Depth</li> </ul>  | : 46 cr  | n       |
| <ul> <li>Height</li> </ul> | : 100 c  | cm      |
| <ul> <li>Weight</li> </ul> | : appro  | x 47 kg |
| <ul> <li>Volume</li> </ul> | : 0,52 r | n³      |
|                            |          |         |





## **GENERAL SPECIFICATIONS**

- Easy cleanable and resistant ABS plastics top platform.
- At both sides there are one piece of quartet and two each quinary drug boxes produced from transparent flexiglass.
- Sliding shelf and drawer for drug preparation.
- There is plastic dustbin one side of trolley and other side serum containers made of ABS plastic.
- Height adjustable IV stand.
- Epoxy painted steel construction is, resistant to corrosion.
- Moves easily by 100 mm diameter four castors two of them with brake to fix the position.

#### **TECHNICAL SPECIFICATIONS**

| <ul> <li>Width</li> </ul>  | : | 93 cm               |
|----------------------------|---|---------------------|
| <ul> <li>Depth</li> </ul>  | : | 46 cm               |
| <ul> <li>Height</li> </ul> | : | 100 cm              |
| <ul> <li>Weight</li> </ul> | : | approx 47 kg        |
| <ul> <li>Volume</li> </ul> | : | 0,52 m <sup>3</sup> |