TROLLEY EMERGENCY AND DRUG DISTRIBUTION 90108302



GENERAL SPECIFICATIONS

- Easy cleanable, strong ABS plastic upper platform and shelf for preparing medicine.
- Useful four pieces drawers can be locked.
- Stainless steel stand defibrillator.
- Stainless steel body.
- Long lived rail system.
- Easy moveable with 125 mm diameter casters and fixed each position with two casters having brake mechanisms.

TECHNICAL SPECIFICATIONS

Width	:	74,5 cm
Depth	:	63 cm
Height	:	105 cm
Weight	:	approx 77 kg
Volume	:	0,61 m ³

