

SH100

Portable Ventilator

ETERNITY



Technical specification

Physical Specification

Dimension: 58mm (H) * 28mm (W) *41 mm (D)

Weight: 15kg (Gross Weight)

Main Unit

Pneumatic Driven Pneumatic Control

Applicable for Adult, Pediatric.

Ventilation: IPPV, Manual(optional) , Pressure Trigger (optional)

Alarm

Low gas source pressure alarm.

Standard Configuration

2L Cylinder

Expiratory valve

1L test bag

Mask

Saddle bag

Silicon tube

Specification

Minute volume	3~20L/min, continuously adjustable	Safe pressure max	≤6Kpa
Breath frequency	10~35bpm, continuously adjustable	Noise	≤65dB (A)
I: E	1: 1.5	Cylinder	2L
O2 concentration	100% / 60%	Gas source pressure:	0.3~0.6Mpa (Flow≥ 50L/min)

BEIJING ETERNITY ELECTRONIC TECHNOLOGY CO., LTD

Add : No .17 Xijing Road , Badachu , Zhonguancun High-tech Zone

Shijingshan District , Beijing 100041 P.R China

TEL: 0086 - 10 - 68827491

FAX : 0086 - 10 - 68817023

www.ernitymed.com

The logo for ETERNITY features the word "ETERNITY" in a bold, blue, sans-serif font. A red, stylized crescent shape is positioned above the letter "N".