# 90101137SS VIGOROUS 11SS HOSPITAL BED 4 MOTORS

## **FUNCTIONS**



Backrest movement by motor 0 to 68°



Kneerest movement by motor 0 to 27 °



Foot rest mechanical adjustment : 0 - 22°



Vascular position



Fowler position



Cardiac position



Trendelenburg position: 0 - 12°



Anti - trendelenburg position: 0 - 12°



Mattress platform minimum height 42 cm



Mattress platform maximum height: 77 cm

## **OPTIONS**

GROUP 1 / CASTOR OPTIONS









**GROUP 4 / BACKREST MOTOR OPTIONS** 

GROUP 2 / SIDE RAIL OPTIONS





GROUP 5 / LEGREST MOTOR OPTIONS





**GROUP 7 / MATTRESS PLATFORM OPTIONS** 





A STATE OF

**GROUP 8 / OTHER ACCESSORY OPTIONS** 

**GROUP 6 / CONTROL PANEL OPTIONS** 



00000





10 40

C3





**DESCRIPTION OF OPTION CODES** 

- TWIN CASTOR TWO WITH BRAKE (120 MM)
- CASTOR CENTRAL & DIRECTIONAL LOCK (125 MM)
- CASTOR CENTRAL & DIRECTIONAL LOCK (150 MM)
- WITHOUT SIDE RAIL
- K5 WITH PE SIDE RAIL
- WITHOUT ELECTRICAL ACCESSORIES
- BATTERY RECHARGEABLE LEAD GEL TYPE. OUTPUT VOLTAGE: 24 V DC. CAPACITY: 1,2 Ah.
- NURSE CONTROL PANEL. PROTECTION CATEGORY : IP 66. 15 FUNCTION KEY. E3 UNDER BED LIGHTING SYSTEM
- FOOT CONTROL PEDAL. ONE PEDAL CPR, ONE PEDAL EXAMINATION POSITION FROM EACH SIDE.
- BACKREST MOTOR (4000 N)
- BACKREST CPR MANUAL CONTROL MOTOR M2
- LEGREST MOTOR (4000 N)
- LEGREST CPR MANUAL CONTROL MOTOR (2000 N)
- Must select only one of them.
- Must select minumum one. Can select 1-2-3 ... both.

- AND REMOTE CONTROL. ROTECTION CATEGORY : IP 66
- REMOVABLE 4 PARTED MESH PLATFORM.
- REMOVABLE DOUBLE WALLED PE BED PLATFORM
- WITHOUT OTHER ACCESSORIES
- IV POLE: DOUBLE HOOK STAINLESS STEEL IV POLE.
- LIFTING POLE : BEARING CAPACITY 75 Kg
- (-RAY AVAILABLE FOR BACK SIDE WITH CASSETTE HOLDER AND COMPACT LAMINATE
- TELEVISION STAND FOR 17" LCD TV
- TRACTION APPARATUS
- BED EXTENSION: 30 cm.
- RESTRAINT STRAP
- W12 CHART HOLDER

# Do not select the others if you select without accessory. If you select the others you can select all or 1-2-3 ... of them

# **GENERAL SPECIFICATIONS**

- General hospital use medical and electrical bed.
- Sections:

Totaly 4 section as below.

- 1-Back rest.
- 2-Seat rest
- 3-Thigh rest.
- 4-Foot rest.
- Design:
- 1-2-3 mm metal tube heavy duty metal construction
- 2- Without space between mattress and head board.
- 3- Stabil low bed scissor construction.
- 4- Epoxy panited metal construction after nano technology
- surface treatment which is easy to clean with disinfectants.
- 5- German Dewert brand motors.

#### Mattress Support type:

- 1-Back rest heavy duty steel grid, optional compact laminate with xray cassette holder.
- 2-Seat rest heavy duty steel grid.
- 3-Thigh rest heavy duty steel grid.
- 4-Foot rest heavy duty steel grid.

#### ■ Movements:

- 1-Backrest movement by motor 0 to 68°.
- 2-Kneerest movement by motor 0 to 27°
- 3-Height adjustment by motor from 42 cm to 77 cm.
- 4-Trendelenburg from 0 to 12°.
- 5-Anti trendelenburg from 0 to 12°.

#### ■ Side rails:

4 parted plastic two pairs for each side.

Strong heavy duty design.

Easy to operate up and down.

C2 right side siderails in and outside dual control panel.

C3 left side, siderails in and outside dual control panel.

Side rails control are easy to clean and disinfect flat surface and strong material. The side rails all spaces design according to European patient safety standart as

EN 60601-2-52. The side rails cover the sleeping platform to prevent patient falling according to

EN 60601-2-52. The side rails locks are safe when it is locked and testes according to EN 60601-2-52.

#### Overall Dimension:

External from bumpers: 102 x 221 cm.

Height adjustment from 42 cm to 77 cm.

# ■ Castors:

150 mm central braking, antistatic and one of them directional lock medical Taiwans Haion brand castors.

# ■ Bumpers:

4 soft toch corner bumpers to protect the bed and environments.

# ■ Head and foot boards:

Plastic heavy duty blow molded easy removable panels. 20 different laminate color options at page 84.

#### ■ Iv pole:

Double hook stainless steel Iv pole.

#### ■ Safe working load:

250 kg at any position.

#### ■ CPR:

With 2 clear visiable red handles on both sides. Writes clearly Cpr on the both handles.

### ■ Bed exstension:

30 cm.

# ■ Drainage bag hooks:

Each longiditual side of the bed.

#### ■ Power supply:

220V / 50 Hz.

With battary backup.







# **DOLSAN**

