

Pregaid


Wireless Basal Thermometer



wearcare


Product feature

Easy to use
High precision 0.05°C
Intelligent analysis
Period prediction&reminder
Ovulation prediction&reminder
Increase pregnancy chance
Charting
Customized reminder
Safe materials
Cloud storage
Partner sharing



Turn on the app

Start the intelligent ovulation prediction



Be aware of body status

PregAid could collect and record users' sublingual temperature from bluetooth smart thermometer everyday, which would be of great diagnostic value.

In-time menstrual period reminder and prediction

PregAid would track user's temperature variations and could predict and remind the user's ovulation and menstrual period correctly.

Choose the best time to conceive

The considerate PregAid would analyze the highest pregnancy possibility and remind the user to get ready to have baby, when the egg is most fresh and mature to create a healthy baby.

Parameter

Dimension	115.6mm(L)×27.1mm(W)×8.3mm(H)
Weight	6.25g(with battery)
Range	25.00°C ~ 45.00°C
Precision	±0.05°C (35.00°C ~ 38.50°C);±0.1°C (25.00°C ~ 34.99°C,38.51°C ~ 45.00°C)
Unit	°C" or"°F" (optional)
App Platform	
IOS	Above iPhone 4s, iPhone, iPhone5s, iPad3, iPad mini
Android (above 4.3)	Cell phone, tablet PC
Battery	
Type	CR1220, 3V/38mAH button cell
Working hours	≥3months (8h/day)

BIOLIGHT CO.,LTD.

Dir: innovation first road, technology innovation coast zhuhai, CHINA Postcode: 519085
Tel: +86-756-3399999 3399983 Fax: +86-756-3399911 Email: overseas@blt.com.cn

 **BLT BIOLIGHT**
<http://www.blt-com.cn>