



# YDW-II

FAR INFRARED HEATER

Stock code [300314]





## Brief Introduction

Ceramic plate as the far infrared radiation resource, and its output can be adjusted;  
Angle of radiant box can be adjusted;  
Height of radiant box can be adjusted;  
Wheels can be locked.

## Standard Configuration

Radiant box, knob with adjustable angle, lift pole, wheels.

## Specification

Power requirement: AC220V-230V/50Hz  
or AC110-120V/50-60Hz, 1000VA

Infrared wave length: 5um~10um

adjustable of heating power: 500W/1000W two levels

Speed of temperature rising:

Under half power:  $\geq 0.18^{\circ}\text{C}/\text{min}$

Under full power:  $\geq 0.36^{\circ}\text{C}/\text{min}$

Angle of radiant box can be adjusted

Adjustable height of radiant box: 1300mm~1700mm



## Environmental Requirement

Environment temperature:  $10^{\circ}\text{C}\sim 20^{\circ}\text{C}$

Wind speed of environment:  $< 0.3\text{m/s}$

## Product Packaging

Package: The whole unit is packed into one carton

Dimension: 955mm x 655mm x 250mm

Gross weight: 25kg



Ningbo David Medical Device Co.,Ltd

Add: Bldg No. 5 Ningbo Smart Park, #98 Chuangyuan Rd, Hi-tech Zone, Ningbo, China

Tel: 0086-574-8780 0008 or 8780 0007

/Fax: 0086-574-8780 3714

E-mail: sales@nbdavid.com

/Web.: www.nbdavid.com

/P.C.: 315048