

ZONCARE-i50

Full Digital B/W Ultrasound Machine

Outstanding performance
Innovation experience



WUHAN ZONCARE BIO-MEDICAL ELECTRONICS CO.,LTD.

ADD: ZONCARE BUILDING, NO. 380, HIGH-TECH 2ND ROAD, EASTLAKE HIGH-TECH
DEVELOPMENT ZONE, WUHAN, HUBEI 430206 P. R. CHINA.

TEL: +86-27-8777 0581 FAX: +86-27-8777 0203

EMAIL: info@zoncare.com <http://www.zoncare.com>

© 2017 WUHAN ZONCARE BIO-MEDICAL ELECTRONICS CO.,LTD. ALL RIGHTS RESERVED. SPECIFICATIONS SUBJECT TO CHANGES WITHOUT
PRIOR NOTICE.



<http://www.zoncare.com>



ZONCARE-i50

Full Digital B/W Ultrasound Machine

- 12.1 inch high resolution medical monitor with 90 viewing angle
- Exclusive employ double wave form imaging technology
- Outstanding portable model, compact and streamlined design
- Double activated transducer connectors, transducer auto-freeze protection
- All-round adjusting knob, efficient and fast
- One key full screen design, hidden menu item, let your diagnosis more intuitive
- Support multiple measurement software package
- Built-in lithium battery, optional trolley and suitcase
- Support multiple data transmission interface
- Auto-generation BMP format report, support USB import and export. Non use any special software to read or print on PC
- Optional laser printer can print the generated report

Outstanding performance, Innovation experience




Trolley




Suitcase

Convenient workflow, excellent imaging quality



High density and broadband frequency transducers



Convex

Linear

Transvaginal

Micro-convex